


KAPA HyperExome V2 Probes

High-performance capture of targets from the most relevant genomic databases.



Efficient whole-exome sequencing (WES) begins with expertly designed probes that effectively capture challenging genomic regions. KAPA HyperExome V2 Probes, the latest Roche WES design, delivers superior coverage of recent versions of ACMGv3.1, RefSeq, CCDS, ClinVar, Ensembl, and COSMIC genomic databases. This design also reduces sequencing requirements, with a compact capture target of 43.2 Mb covering 37.5 Mb of key content.

- **Leverage the most relevant genomic databases**
- **Unlock unique insights in clinical research** by capturing more content from key genomic databases
- **Achieve higher uniformity** across the entire range of %GC content
- **Sequence with confidence** with an improved whole-exome workflow
- **Streamline targeted sequencing** with the HyperCap Workflow, driven by KAPA HyperPrep, KAPA HyperPlus, and KAPA EvoPlus Library Prep Kits



Leverage the most relevant databases

- Build confidence by using a design featuring updated ClinVar, RefSeq & Ensembl content*
- Achieve exceptional coverage of ACMGv3.1, CCDS, ClinVar, COSMIC, Ensembl, HGMD, & RefSeq databases
- Balance specificity & coverage with probes designed using the telomere-to-telomere genome assembly CHM13




Figure 1. Superior probe database coverage by the KAPA HyperExome V2 design compared to other suppliers' designs across important genomic databases (better coverage by as much as 17%). Database data retrieved in Jan 2023, unpadded capture target used to compare across suppliers. The KAPA HyperExome V2 panel was designed to cover coding exon sequence from the following annotation sources: RefSeq (June 29, 2022), ClinVar (June 29, 2022) and Ensembl release v106.

Unlock unique insights in clinical research

- Cover more bases and reduce “blind” spots, leaving fewer unknowns
- Increase value and confidence in the results with superior coverage of database content




Figure 2. Percent bases covered by at least 30X across some key genomic databases. KAPA HyperExome V2 Probes deliver better database coverage compared to Suppliers I, A, and T, which leave more “blind” spots and may thus require repeated experiments and/or the addition of supplemental probes. Supplier sample prep protocols were followed with singleplex 4h hybridizations for Supplier I (3 replicates, NA12878), 8-plex O/N hybridizations for Supplier A (48 data points from 6 replicate captures of 16 coriell DNAs), 8-plex O/N hybridizations for Supplier T (72 data points from 16 coriell and 24 blood-extracted DNAs) and 16-plex O/N hybridizations for the KAPA HyperExome V2 Probes (48 data points from 16 blood-extracted DNAs in 3 captures) with the KAPA EvoPlus Kit, KAPA HyperPlus Kit, and KAPA HyperPrep Kit (48 data points from 16 coriell DNAs in 3 captures). KAPA HyperExome V2 Probes-enriched libraries were sequenced on a NovaSeq™ 6000 System at 2 x 100 bp and 60 M high-quality reads were analyzed per library. Supplier-enriched libraries were sequenced proportionally to their capture target size.

*(end June 2022)

For Research Use Only. Not for use in diagnostic procedures.

Achieve higher uniformity across the entire %GC range

- Obtain exceptionally uniform coverage even through the extremes of the %GC spectrum
- Eliminate GC bias to cover low- and high-GC regions equally well, with an optimized design and uniform library amplification




Figure 3. Exceptional uniformity of normalized capture coverage across the extremes of the GC% spectrum. The blue dashed line represents the optimal uniformity in the ideal state in which all regions—regardless of their GC content—would be equally covered. Exome libraries were prepared using KAPA HyperExome V2 Probes, Supplier T exome probes, or Supplier A exome probes as described in Figure 2. Sequencing and analysis were carried out as described for Figure 2.

Sequence with confidence


Increase result confidence by following sample identities throughout the workflow with sample-tracking endogenous SNPs covered by the KAPA HyperExome V2 Probes. Or, explore new potential capabilities with a unique set of 96 probes that are included in the panel and designed to capture non-naturally occurring sequences that can be spiked into samples.

- Select from an extensive list of 529 sample-tracking SNPs that includes the Pengelly¹ and Yousefi² sets
- Rely on high precision (99.49%) and recall (98.95%) for SNP detection
- Explore new possibilities with 96 utility probes embedded in the design, which may be used to capture exogenous synthetic DNA fragments, potentially used as process controls

Streamline targeted sequencing with the KAPA HyperCap Workflow

KAPA HyperCap Workflow delivers complex libraries by combining the high conversion rate of KAPA HyperPrep, KAPA HyperPlus, and KAPA EvoPlus Kits with KAPA HyperCap probes (such as KAPA HyperExome V2 Probes), creating a streamlined, single-vendor supported workflow.

- Multiplex up to 16 samples in the same capture, and potentially post-capture multiplex more samples in the same sequencing lane, with KAPA UDI Primer Mixes, 1-384
- Reduce workflow complexity and hands-on time with KAPA Universal Enhancing Oligos, eliminating the need for adapter-matched blocking oligos
- Automate the entire KAPA HyperCap Workflow without the need for a SpeedVac—now with all hybridization and bead wash steps at 55°C



1. A SNP profiling panel for sample tracking in whole-exome sequencing studies. Pengelly RJ, et al. *Genome Med.* 2013 Sep 27;5(9):89. doi: 10.1186/gm492. eCollection 2013

2. A SNP panel for identification of DNA and RNA specimens. Yousefi S, et al. *BMC Genomics.* 2018 Jan 25;19(1):90. doi: 10.1186/s12864-018-4482-7.

Project: Roche KAPA HyperExome V2 Launch Data, Pleasanton, CA, Mar. 2023

For Research Use Only. Not for use in diagnostic procedures.

Ordering Information

Product name	Capture reactions	Catalog No.
KAPA HyperExome V2 Prep Kit Complete KAPA HyperCap Workflow NGS sample prep kit for 192 8-plex reaction samples <i>Includes: KAPA HyperPrep Kit, Universal Adapters, HyperBeads, Probes Resuspension Buffer, HyperExome V2 Probes, HyperCapture Reagent Kit, and Hypercapture Bead Kit</i>	192 8-plex rxn	9983759001
KAPA HyperExome V2 Plus Kit Complete KAPA HyperCap Workflow NGS sample prep kit for 192 8-plex reaction samples <i>Includes: KAPA HyperPlus Kit, Universal Adapters, HyperBeads, Probes Resuspension Buffer, HyperExome V2 Probes, HyperCapture Reagent Kit, and Hypercapture Bead Kit</i>	192 8-plex rxn	9983775001
KAPA HyperExome V2 Evo Kit Complete KAPA HyperCap Workflow NGS sample prep kit for 192 8-plex reaction samples <i>Includes: KAPA EvoPlus Kit, Universal Adapters, HyperBeads, Probes Resuspension Buffer, HyperExome V2 Probes, HyperCapture Reagent Kit, and Hypercapture Bead Kit</i>	192 8-plex rxn	9983783001
	12	9718630001
	24	9718648001
	48	9718656001
	96	9718664001
KAPA HyperExome V2 Probes	192	9718672001
	384	9718699001
	768	9718702001
	1152	9718729001
	1536	9718737001

Demo data

Demo data files are available for evaluation. Please contact a local Roche representative.

KAPA HyperExome V2 Probes

Learn more at go.roche.com/HyperExomeV2

KAPA HyperCap Workflow

Learn more at go.roche.com/HyperCap

Published by:
Roche Sequencing and Life Science
9115 Hague Road
Indianapolis, IN 46256

sequencing.roche.com

AVENIO, KAPA, KAPA EVOPLUS, KAPA HYPERCAP, KAPA HYPEREXOME, KAPA HYPERPLUS, and KAPA HYPERPREP are trademarks of Roche. All other product names and trademarks are the property of their respective owners.
For Research Use Only. Not for use in diagnostic procedures.
© 2023 Roche Sequencing and Life Science. All rights reserved.

