

Inter-site workflow performance of a fully automated sample preparation system for NGS applications: *integration with KAPA HyperExome probes* for whole exome sequencing, the AVENIO Edge† instrument, and navify® Mutation Caller#*

Rachel Kasinskas¹, Yanpeng Xi², Lara Sucheston-Campbell³, Pierre-Luc Janvier², Helle Sorensen³, and Alexandra Markus²

¹Roche Sequencing and Life Science, Wilmington, MA

²Roche Diagnostics Solutions, Pleasanton, CA

³Roche Diagnostics Solutions, Santa Clara, CA

³ Lara Sucheston-Campbell is no longer an employee of Roche.

navify Mutation Caller is not for sale/use in the United States. † The AVENIO Edge System and navify® Mutation Caller are for Research Use Only workflows. * Reagents and kits mentioned are for Research Use Only. Not for use in diagnostic procedures.

MC-US-12966

Background and Goals

- Next generation sequencing (NGS) sample preparation automation systems, for library preparation and target enrichment, are used by laboratories to drive scalability in precision medicine research.¹⁻³ The AVENIO Edge System[†] is a fully automated NGS sample preparation system for library preparation, target enrichment, pooling, normalization and quantification.

This study evaluated the inter-site performance of the AVENIO Edge System[†] automated workflow for whole-exome sequencing (WES) with the KAPA HyperExome Probes* in comparison to the manual workflow at four sites.

1. Results – Part 1

- For the baseline runs, percent reads on-target +/- SD were 87.3±0.5 (A) and 87.4±1.1 (B), mean depth of coverage +/- SD 55.7±1.4 (A) and 54.8±1.7 (B), and percent total duplication +/- SD 4.8±1.5 (A) and 4.9±1.6 (B). Percent reads on-target +/- SD were 86.4±0.7 (A) and 86.3±1.2 (B) [50ng], 87.2±0.4 (A) and 87.2±0.6 (B) [8-plex], and 81.4±1.7 (A) and 83.2±0.9 (B) [manual]. Six manual run (A) samples did not meet acceptance criteria (>60M reads for panel ≥40 Mb).




Figure 1: AVENIO Edge[†] workflow performance and manual workflow comparison: Percent reads on-target (A), mean depth of coverage (B), and percent total duplication (C).

2. Results – Part 2

- Automated workflow recall, precision and concordance was greater than 99.1% (A) and 99.1% (B) (Table 1). Eight samples [six Run6 (A), one Run2 (B), one Run5 (B)] that did not meet acceptance criteria were removed from variant calling analysis since low coverage or base quality can result in wrong calls or a refusal to call.

Table 1: Workflow recall, precision and concordance by navify® Mutation Caller#

		AVENIO Edge [†] Automation					Manual
		Run 1 (baseline)	Run 2 (baseline)	Run 3 (baseline)	Run 4	Run 5	Run 6
Site A	% Recall	99.19	99.25	99.25	99.27	99.28	99.43
	% Precision	99.16	99.21	99.24	99.31	99.24	99.49
	% Concordance	99.10	99.16	99.19	99.26	99.19	99.45
Site B	% Recall	99.24	99.15	99.23	99.26	99.26	99.36
	% Precision	99.20	99.14	99.22	99.27	99.30	99.35
	% Concordance	99.13	99.07	99.16	99.22	99.24	99.30

Acknowledgements

Bill LaRochelle, Persis Wadia, Lisha Capucion, Alva Liang, Shannon Rich, Malek Kunbargi, Michael Pintor, Marsha McMakin, Jenny Gu, Mauricio Salvo, Markos Mihalatos, Nripesh Prasad, Shawn Levy, Bernd Timmermann

Disclosures

Sponsored by Roche Sequencing Solutions (Pleasanton, CA, USA), part of the Roche organization, in collaboration with Discovery Life Sciences, Huntsville, AL, USA and Max Planck Institute, Berlin, Germany. AVENIO EDGE, KAPA HYPERCAP, KAPA HYPEREXOME and NAVIFY are trademarks of Roche. All other product names and trademarks are the property of their respective owners.


Abbreviations

gDNA, genomic DNA; NA, nucleic acid; Rep, replicate; SD, standard deviation; SNP, single nucleotide polymorphism; WES, whole exome sequencing

3. Results – Part 3

- Consistent inter-instrument reproducibility was observed across each instrument site (A/B and two Roche sites; seven runs, five instruments) (Figure 2). All AVENIO Edge[†] instruments achieved >85% reads on-target (Fig 2A & 2C). Individual read distribution for 24 replicates across 7 instrument runs (168 samples) showed average percent reads on-target ± SD, 88±0.57 (A) (Fig 2B) and 86±1.62(B) (Fig 2D).

Figure 2: Inter-site and inter-instrument reproducibility: Roche Pleasanton (3 AVENIO Edges[†], 3 runs), Roche Wilmington (1 AVENIO Edge[†], 1 run) and independent labs [1 AVENIO Edge[†] per lab, 3 runs(1-3)], A (A & B) and B (C & D).



Methods

- Two independent labs (site A and B) each prepared 24 replicates (NA12878 human gDNA) for each run (Table 2). Five automated runs evaluating nucleic acid inputs (100 ng and 50 ng) and pooling (singleplex and 8-plex) conditions were conducted following AVENIO Edge System[†] instructions for use with KAPA HyperExome probes* (Roche) and compared against a (manual) KAPA HyperCap v3.0 workflow* (Roche) run (100 ng, singleplex) with Illumina sequencing. At each lab, repeatability and well-to-well comparison was measured across three baseline instrument runs (100 ng, singleplex) using percent reads on-target, mean depth of coverage, and percent total duplication and standard deviation (SD) of each measure. Percent reads on-target were also assessed for the 50ng and 8-plex pooling instrument runs, and the manual run. Variant calling was performed with navify® Mutation Caller# (Roche) to determine the variant detection performance of the automated workflow as compared to NA12878 gold standard with 24,901 high confidence SNPs in the KAPA HyperExome* target region. Inter-site and inter-instrument reproducibility was evaluated against two Roche sites.

Table 2. Study design: AVENIO Edge[†] and manual study runs and conditions tested by each independent lab

Run	Workflow	Sample Type	Replicates	Nucleic Acid Input	Pooling	Description
1	Automated	NA12878	24	100 ng	Singleplex	Rep1 (baseline)
2			24	100 ng	Singleplex	Rep2 (baseline)
3			24	100 ng	Singleplex	Rep3 (baseline)
4			24	100 ng	8-plex	Pre-capture pooling performance
5			24	50 ng	Singleplex	Low NA sample input performance
6	Manual		24	100 ng	Singleplex	Compare AVENIO Edge [†] automated runs to a manual run

Conclusions

- AVENIO Edge System[†] for NGS sample preparation automation integrated with navify® Mutation Caller# for variant calling demonstrated robust inter-site workflow and variant detection performance as compared to a manual workflow for WES using KAPA HyperExome Probes* to enable improved laboratory efficiency and ensure result confidence.

References

- Collins FS, Varmus H. A new initiative on precision medicine. N Engl J Med 2015; 372:793- 5.
- May M. Automated sample preparation. Science. 2016. 351(6270), 300-302.
- Kong F, Yuan L, Zheng YF et al. Automatic Liquid Handling for Life Science: A Critical Review of the Current State of the Art. J Lab Auto 2012; 17:169-185.