

INSTRUMENT COMPATIBILITY TABLE



- Preferred kit for the instrument
- Can be used, but requires the addition of ROX reference dye
- * Addition of ROX High is required
- † Addition of ROX Low is required

KAPA
Cat. No. Roche
Cat. No.

ABI Prism® 7000
ABI Prism® 7700
ABI Prism® 7900HT
ABI 5700
ABI 7300
ABI StepOne™/StepOnePlus™
ABI 7500
ABI ViiA™ 7
ABI QuantStudio™ Family
Stratagene Mx3000P™
Stratagene Mx3005P™
Stratagene Mx4000™
Bio-Rad iCycler iQ®/iQ® 5
Bio-Rad iCycler MyiQ™
Bio-Rad CFX96™/CFX384™
Bio-Rad Opticon™
Bio-Rad Chromo4™
Cepheid SmartCycler®
Corbett Rotor-Gene™
Qiagen Rotor-Gene™ Q
Eppendorf Mastercycler®
Illumina Eco™
Roche LightCycler® 480/1536
Roche LightCycler® 96

Universal – Includes separate tubes of ROX passive reference dye (High and Low)

Kit	KAPA Cat. No.	Roche Cat. No.	Volume	ABI Prism® 7000	ABI Prism® 7700	ABI Prism® 7900HT	ABI 5700	ABI 7300	ABI StepOne™/StepOnePlus™	ABI 7500	ABI ViiA™ 7	ABI QuantStudio™ Family	Stratagene Mx3000P™	Stratagene Mx3005P™	Stratagene Mx4000™	Bio-Rad iCycler iQ®/iQ® 5	Bio-Rad iCycler MyiQ™	Bio-Rad CFX96™/CFX384™	Bio-Rad Opticon™	Bio-Rad Chromo4™	Cepheid SmartCycler®	Corbett Rotor-Gene™	Qiagen Rotor-Gene™ Q	Eppendorf Mastercycler®	Illumina Eco™	Roche LightCycler® 480/1536	Roche LightCycler® 96
KAPA SYBR® FAST qPCR Kit	KK4600	07959362001	1 x 1 mL	○*	○*	○*	○*	○*	○*	○†	○†	○†	○†	○†	○†			●	●	●	●	●	●	●	●	●	●
	KK4601	07959389001	1 x 5 mL	○*	○*	○*	○*	○*	○*	○†	○†	○†	○†	○†	○†			●	●	●	●	●	●	●	●	●	●
	KK4602	07959397001	2 x 5 mL	○*	○*	○*	○*	○*	○*	○†	○†	○†	○†	○†	○†			●	●	●	●	●	●	●	●	●	●
KAPA SYBR® FAST One-Step qRT-PCR Kit	KK4650	07959613001	1 x 1 mL	○*	○*	○*	○*	○*	○*	○†	○†	○†	○†	○†	○†			●	●	●	●	●	●	●	●	●	●
	KK4651	07959621001	1 x 5 mL	○*	○*	○*	○*	○*	○*	○†	○†	○†	○†	○†	○†			●	●	●	●	●	●	●	●	●	●
	KK4652	07959630001	2 x 5 mL	○*	○*	○*	○*	○*	○*	○†	○†	○†	○†	○†	○†			●	●	●	●	●	●	●	●	●	●
KAPA PROBE FAST qPCR Kit	KK4701	07959800001	1 x 1 mL	○*	○*	○*	○*	○*	○*	○*	○*	○*	○*	○*	○*			●	●	●	●	●	●	●	●	●	●
	KK4702	07959818001	1 x 5 mL	○*	○*	○*	○*	○*	○*	○*	○*	○*	○*	○*	○*			●	●	●	●	●	●	●	●	●	●
	KK4703	07959826001	2 x 5 mL	○*	○*	○*	○*	○*	○*	○*	○*	○*	○*	○*	○*			●	●	●	●	●	●	●	●	●	●
KAPA PROBE FORCE qPCR Kit	KK4300	07959320001	100 x 20	●	●	●	●	●	●	●	●	●	●	●	●			●	●	●	●	●	●	●	●	●	●
	KK4301	07959338001	500 x 20	●	●	●	●	●	●	●	●	●	●	●	●			●	●	●	●	●	●	●	●	●	●
	KK4302	07959346001	1000 x 20	●	●	●	●	●	●	●	●	●	●	●	●			●	●	●	●	●	●	●	●	●	●
	KK4303	07959354001	5000 x 20	●	●	●	●	●	●	●	●	●	●	●	●			●	●	●	●	●	●	●	●	●	●

ABI Prism® – contains ROX High passive reference dye

Kit	KAPA Cat. No.	Roche Cat. No.	Volume	ABI Prism® 7000	ABI Prism® 7700	ABI Prism® 7900HT	ABI 5700	ABI 7300	ABI StepOne™/StepOnePlus™	ABI 7500	ABI ViiA™ 7	ABI QuantStudio™ Family	Stratagene Mx3000P™	Stratagene Mx3005P™	Stratagene Mx4000™	Bio-Rad iCycler iQ®/iQ® 5	Bio-Rad iCycler MyiQ™	Bio-Rad CFX96™/CFX384™	Bio-Rad Opticon™	Bio-Rad Chromo4™	Cepheid SmartCycler®	Corbett Rotor-Gene™	Qiagen Rotor-Gene™ Q	Eppendorf Mastercycler®	Illumina Eco™	Roche LightCycler® 480/1536	Roche LightCycler® 96
KAPA SYBR® FAST qPCR Kit	KK4603	07959419001	1 x 1 mL	●	●	●	●	●	●	●																	
	KK4604	07959427001	1 x 5 mL	●	●	●	●	●	●	●																	
	KK4605	07959435001	2 x 5 mL	●	●	●	●	●	●	●																	
KAPA SYBR® FAST One-Step qRT-PCR Kit	KK4660	07959656001	1 x 1 mL	●	●	●	●	●	●	●																	
	KK4661	07959664001	1 x 5 mL	●	●	●	●	●	●	●																	
	KK4662	07959672001	2 x 5 mL	●	●	●	●	●	●	●																	
KAPA PROBE FAST qPCR Kit	KK4705	07959842001	1 x 1 mL	●	●	●	●	●	●	●																	
	KK4706	07959869001	1 x 5 mL	●	●	●	●	●	●	●																	
	KK4707	07959877001	2 x 5 mL	●	●	●	●	●	●	●																	

ROX Low – contains ROX Low passive reference dye

Kit	KAPA Cat. No.	Roche Cat. No.	Volume	ABI Prism® 7000	ABI Prism® 7700	ABI Prism® 7900HT	ABI 5700	ABI 7300	ABI StepOne™/StepOnePlus™	ABI 7500	ABI ViiA™ 7	ABI QuantStudio™ Family	Stratagene Mx3000P™	Stratagene Mx3005P™	Stratagene Mx4000™	Bio-Rad iCycler iQ®/iQ® 5	Bio-Rad iCycler MyiQ™	Bio-Rad CFX96™/CFX384™	Bio-Rad Opticon™	Bio-Rad Chromo4™	Cepheid SmartCycler®	Corbett Rotor-Gene™	Qiagen Rotor-Gene™ Q	Eppendorf Mastercycler®	Illumina Eco™	Roche LightCycler® 480/1536	Roche LightCycler® 96
KAPA SYBR® FAST qPCR Kit	KK4619	07959575001	1 x 1 mL							●	●	●	●	●	●												
	KK4620	07959583001	1 x 5 mL							●	●	●	●	●	●												
	KK4621	07959591001	2 x 5 mL							●	●	●	●	●	●												
KAPA PROBE FAST qPCR Kit	KK4716	07959940001	1 x 1 mL							●	●	●	●	●	●												
	KK4717	07959958001	1 x 5 mL							●	●	●	●	●	●												
	KK4718	07959966001	2 x 5 mL							●	●	●	●	●	●												

